

Product Usage Guide

All Purpose Sanitiser Unfragranced


Health & Safety Please refer to relevant MSDS

An all-purpose sanitiser based on natural fermentation extracts and mild, biodegradable surfactant. Kills bacteria to EN1276. Approved to EN14476 for enveloped viruses, including coronavirus 2019 n-CoV / SARS-CoV-2.






1




Tip bottle and use guide to measure correct dilution in the measuring chamber. Remove cap and pour.



Dilute 10ml of product into a 750ml trigger spray.



Spray surface, leave for the required contact time, (15 seconds for enveloped viruses, 30 seconds for bacteria) and wipe.**



Tip bottle and use guide to measure correct dilution in the measuring chamber. Remove cap and pour.



Dilute 60ml of product into a 5L bucket.

4

Mop surface, leave for the required contact time, (15 seconds for enveloped viruses, 30 seconds for bacteria). No rinsing necessary. Leave to air dry.**



**please ensure surfaces are clean prior to sanitising